


VideoStar Frenzel Goggles, an ideal complement to NeuroCom® Balance Solutions



A New Step in Observing Eye Movement

Natus offers VideoStar Video Frenzel Goggles, which are ideal for observing eye movement and recording video during Canalith Repositioning maneuvers for BPPV therapy. Video recorded eye movements are also used for teaching, patient education and in medical/legal cases.



VideoStar Monocular (D1 140930)

Comfort: featuring an ultra-light goggle, the VideoStar is one of the most comfortable systems on the market

Quality: using a high quality infra-red camera, the VideoStar requires only one infra-red LED to illuminate the eyeball

Flexibility: the camera can be placed on the right or left side with magnets

Natus — Your Balance Experts


- ▶ Balance technology & software solutions
- ◆ Balance diagnostics
- ◆ Functional assessment & rehabilitation
- ▶ Concussion management & fall prevention
- ▶ Resources to build your balance center

What's Included?

- **VideoStar Video Frenzel Goggles**
- **VideoStar Monocular Camera** features a USB 2.0 camera with high resolution and a high sensitivity to give you clear and high quality video
- **VideoStar DiSoft II Software** compatible with Windows® 8, the **VideoStar software was developed in response to customer requests** with complete nystagmus testing, digitally recorded eye movement videos, video reply and comparison
- **Laptop Computer** 15.6" laptop computer with USB 2.0 ports running Windows 8
- **VideoStar USB Cable** includes five meter USB cable

Options

- Binocular Free-field option allows recording of both eyes simultaneously
- Remote control for processes such as start/stop, recording, next test while staying near the patient
- View the patient's eyes on an extended 32" TV monitor



VideoStar Binocular (D1 140932)

Technical specifications

Specifications

- USB 2.0 camera
- Resolution: 1024 x 768 pixels
- Video frame rate: 50 Hz
- Goggle weight: 102 gr (monocular)

Safety Standards

- IEC 60601-1-1 Patient Safety Standard
- IEC 60601-1-2 Electromagnetic compatibility



VideoStar Binocular (D1 140932)



Distributed by:

natus

Natus Medical Incorporated
1501 Industrial Road
San Carlos, CA 94070 USA
1-800-303-0306
+1-650-802-0400
www.natus.com

For more information, please contact your local Natus District Sales Manager at 1-800-303-0306 or visit www.natus.com

Visit our **NERVE Center®** education portal at nervecenter.natus.com

© 2014 Natus Medical Incorporated.